

PPRB™

Patient and Personal Relief Bag


DISPOSABLE UNISEX URINE BAG

Directions for Use:


1. Use the PPRB



2. Seal the Unisex Collar



3. Dispose in Trash Container



alphapointe™
FACTICAL DIVISION

87-46 123rd Street, Richmond Hill, NY 11418
www.alphapointe.org 718.854.7300

CAGE Code: 056F9 PN: 60049



alphapointe™
FACTICAL DIVISION

87-46 123rd Street, Richmond Hill, NY 11418
www.alphapointe.org 718.854.7300

CAGE Code: 056F9 PN: 60049

100 ml

200 ml

300 ml

400 ml

Quantities are approximations only

To Measure Quantity for Patient Use:

MEDICAL PERSONNEL

1. Wait 60 seconds to form gel
2. Place horizontally on flat surface
3. Gently apply pressure from front to back to remove excess bulging
4. Measure contents by determining closest quantity fill line